

*Work Module 1. Liquid Biopsy and Biomarkers*

Experimental Protocol  
**Plasma preparation for liquid biopsy studies**

**Total time for experiment:** 40 minutes

**Reagents:**

4ml or 10ml Vacutainer K2 EDTA tube (Becton Dickinson, ref: # 367864 (4ml) #367525 (10ml))

Tubos Streck Cell-Free DNA BCT CE 10 ml ( ref: # 218997 (100 tubes))

4ml polypropylene cryogenic vial, round bottom, self-standing (e.g. Corning, #430662) or 1ml Eppendorf tubes

10ml and 5ml serological disposable pipettes or p1000 filter tips

15ml polypropylene centrifuge tubes

Freezer storage boxes for 1-4ml cryogenic vials (e.g. Corning, #431120)

**Instruments and tools:**

Centrifuge, capable of ~3000g with a swing bucket rotor

Pipetting aid or plastic disposable bulb pipettes or p1000 pipette

-80°C Temperature Freezer

**Prior considerations:**

Samples should ideally be processed within 4 hours from the time of extraction when using EDTA tubes or 24-48 if using Streck tubes.

Ideally Streck tubes should be used for blood collection as they contain a preservative that stabilizes nucleated blood cells and prevents the release of contaminating genomic DNA. In theory, cfDNA and gDNA are stable for up to 14 days at 6 °C to 37 °C after blood extraction. However, the samples should be processed within 24-48 hours in order to avoid contamination with lysed white blood cells.

**Steps:**

- Record the time and date of the sample extraction as well as the subject ID in an appropriately designed sample database.

- Centrifuge the EDTA tube at room temperature (15°C to 25°C) for 10 min at 1600 ( $\pm$  150)g.

*N.B. sample in Streck tubes can be centrifuged at 3000g for 20 minutes.*

- After centrifugation, carefully remove tube from centrifuge and transfer the plasma (supernatant) to a 15ml centrifuge tube without disturbing the cellular layer using a disposable serological pipette (or disposable bulb pipette or p1000 pipettes with filter tip). *N.B. Aspirate ~5mm above the leukocyte layer to reduce contamination of the plasma with cells.*
- Centrifuge the plasma in the 15ml centrifuge tube at room temperature (15°C to 25°C) for 10min at 3000 ( $\pm$ 150)g to remove any residual intact blood cells carried over from the first centrifugation step.

*N.B. sample can also be centrifuged at 3000g for 20 minutes.*

- After centrifugation, carefully remove tube from centrifuge and transfer 1ml-4ml of plasma to 1-4ml polypropylene cryogenic vials using a disposable serological pipette (or disposable bulb pipette or p1000 pipettes with filter tip). *N.B. Leave a residual volume of approximately 0.3 ml (~7mm) on the bottom of the 15ml tube to avoid contaminating the plasma with cells.*
- Immediately freeze plasma upright in storage box at  $-80^{\circ}\text{C}$
- Collect the cellular layer (buffy coat) using a P1000 and filtered tip and transfer it to a 2 mL Eppendorf. Immediately freeze it and storage in  $-20^{\circ}\text{C}$  or  $-80^{\circ}\text{C}$  freezer.

*N.B. Aliquots of an adequate volume should be used in order to avoid freeze-thawing of the samples. The volume and number of aliquots will depend on the downstream use of the samples and also the storage capacity of the  $-80$  freezer. 1ml aliquots are useful as this is the standard volume used for many cfDNA extraction protocols and various aliquots can be defrosted if more volume is required.*

**Contact:**

liquidbiopsy@ciberonc.es