

Work Module 1. Liquid Biopsy and Biomarkers

Experimental Protocol

Platelets preparation for liquid biopsy studies

Total time for experiment: 40 minutes

Reagents:

4ml or 10ml Vacutainer K2 EDTA tube (Becton Dickinson, ref: # 367864 (4ml) #367525 (10ml))

2 mL Eppendorf tubes

1ml or 3ml serological disposable pipettes or p1000 filter tips

15ml polypropylene centrifuge tubes

Freezer storage boxes for 1-4ml cryogenic vials (e.g. Corning, #431120)

RNAlater™ Stabilization Solution (e.g. Thermo Fisher, #AM7020)

Instruments and tools:

Centrifuge, capable of ~3000g with a swing bucket rotor

Pipetting aid or plastic disposable bulb pipettes or p1000 pipette

-80°C Temperature Freezer

Prior considerations:

Samples should ideally be processed within 4 hours from the time of extraction when using EDTA tubes.

Steps:

- Record the time and date of the sample extraction as well as the subject ID in an appropriately designed sample database.
- First centrifugation of EDTA tubes at low r.p.m (speeds), 120g for 20 minutes at room temperature (RT).

- After this centrifugation, transfer 9/10 of plasma (containing platelets) to the 15 mL falcon tube.
- Centrifuge the Falcon tube with the 9/10 of plasma at 360 g for 20 minutes.
- Transfer supernatant to a new 15 mL falcon tube labeled as plasma without touching the pellet with platelets.
- Resuspend the pellet of platelets with the residual plasma and transfer platelets to the 2 mL tubes with screw caps.

Add 60 to 80 µL of RNAlater solution to the tube and mix gently.

- Refrigerate tubes at 4°C overnight and freeze at -80°C the following day.

N.B. Aliquots of an adequate volume should be used in order to avoid freeze-thawing of the samples. The volume and number of aliquots will depend on the downstream use of the samples and also the storage capacity of the -80 freezer. 1ml aliquots are useful as this is the standard volume used for many cfDNA extraction protocols and various aliquots can be defrosted if more volume is required.

Contact:

liquidbiopsy@ciberonc.es